

Product Code . DS-LE-10905

Top Loading Balance (Digital)



Description

DESCRIPTION :-

Auto shut off.

Auto calibration.

Auto zero tracking.

Auto backlight.

Low battery indication.

Kg/g/oz conversion.

Large LCD (6 digits 15 mm height).

Large Square weighing pan.

Net weight/stability indication.

1/2000, 1/5000 & 1/15000 division available.

Highly efficient

Precise result

Bright LED display

Maximum weight: Min 1200 g

Minimum unit weight: 0.1 g

Linearity: $\pm 0.003\text{g}$

Maximum Load: Min 0.02g

Repeatability/Std. Dev.: Min 0.001g

Sensitivity Drift: $\pm 20\text{ppm}/^{\circ}\text{C}$ ($10^{\circ}\text{C}\sim 30^{\circ}\text{C}/50^{\circ}\text{F}\sim 86^{\circ}\text{F}$)

Display Type: 7 segment LCD with backlight (character height 16mm)

Display Refresh: 10 times/second,

Pan size (W) x (D): 133X170 mm

Display Modes: g & oz

Key function: ON/TARE, MODE, OFF, LCD display

Auto / Manual Shut Off

Standard accessories with operation manual

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/