



Product Code . DS-LE-10906

Distillation Apparatus

Description

Purifying Globally: Glass Distillation Unit Exporter from the UK

The UK is a leading source for high-quality **Glass Distillation Units Exporter UK**, renowned for their precision and reliability. These trusted exporters deliver a range of units suitable for various applications, from purifying laboratory-grade water to producing essential oils and other vital solutions. Explore the UK glass distillation unit export landscape and find the perfect partner to meet your specific purification needs.

DESCRIPTION :-

Stainless steel without water-controlled electric heating distilled water devices to tap water as source of water, electric heating, distillation system, pure water, for health care, pharmaceutical, plating, and laboratory and other departments to use, this product rational structure, easy to use, the operation simple, out of water, large, safe, reliable, stable quality, durable.

The products should all selected highquality stainless-steel material

The use of advanced technology means well made.,

Condenser plate stainless steel seamless pipe system, heat exchange efficiency is high, a large amount of water

Made of high quality stainless steel materials & by advanced technology

Condenser should be made of seamless stainless steel pipes for very good, cooling effects & auto system.



Immersing electric heating tube

Capacity : 10 L/h (Min)

Power: 7.5 kw (Min)

Power Supply: 380 V AC

Standard accessories with operation manual

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/