



Product Code . DS-LE-10908

pH Meter

Description

DESCRIPTION :-

This digital pen style pH meter is an ideal instrument for anywhere you may need to test pH levels.

Very accurate and durable

Professional PH/mV METER

pH range : 0 to 14 PH mV range : 0 - 1800 mV

Wide manual temperature compensation adjustment can be operated by push button on the front panel easily

Calibration for pH 7, pH 4 and pH 10 by push button



Auto shut off saves battery life, Memorize Max/Min reading with recall

Max pH : 0.02 pH

Dimension: min 219 x 199 x 84mm

Large LCD, dual function display

Wide applications: water conditioning, aquariums, beverage

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/