



Product Code . DS-LE-10909

Kipp's Apparatus (Complete Set)

Description

DESCRIPTION :-

Capacity: 250ml

Heat resistant

Chemical Resistant



Usage/Application: Laboratory

Material : Borosilicate Glass

Color : Transparent

With Interchangeable joint

Kipp Hydrogen Gas Generator should a typical device for producing gas by reacting solid and liquid reagents under normal temperature conditions (without heating). Such as hydrogen, carbon dioxide, hydrogen sulfide, etc. should be used to prepare it

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/