

Product Code . DS-LE-10912

Conical Flask with Side Table



Description

DESCRIPTION :-

Color: Transparent

Material: Borosilicate Glass

Capacity: Min 250ml, Height: Min 155mm, Body Dia: Min 85 mm, Neck I.D.: Min 35 mm

Capacity: Min 500ml, Height: Min 185mm, Body Dia: Min 105 mm, Neck I.D.: Min 35 mm

Shape : Cylindrical

Graduated

Heat resistant

Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/