

Product Code . DS-LE-10913

Weighting Bottle



Description

DESCRIPTION :-

Capacity: Min 150ml

Bottle height: Min105mm

Body Dia: Min 62mm

Thread size(mm): Min GL45



Material: Borosilicate Glass

Bottles complete with polypropylene cap and drip-free pouring ring.

Complies with ISO 4796., Ideal for mixing and storage of culture medium.

Sloping shoulders for smooth pouring action.

For use up to 140°C

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/