



Product Code . DS-LE-10914

Funnel with Long Tube

Description

DESCRIPTION :-

Made of best quality & heat proof

Shape: Conical

Material: Borosilicate glass

With 60° bowl angle for precise fitting of filter paper



Diameter: Min 90 mm

Steam O.D: Min 11mm

Steam length: Min 90mm

Heat resistant

Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/