



Product Code . DS-LE-10915

Porcelain Basin

Description

DESCRIPTION :-

Material : Porcelain

Color : White

Glazed finish

Chemical resistant

Round bottom with spout & glazed

Capacity 60 ml

? Outside dia : 70 mm

? Height : 28 mm

Capacity : 75 ml,

? Outside dia : 80 mm

? Height : 30 mm ,

Capacity : 100 ml,

? Outside dia : 88 mm

? Height : 34 mm ,

Capacity : 150 ml,

? Outside dia : 100 mm ,

? Height : 40 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/