

Product Code . DS-LE-10916

Test Tube



Description

DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass

With rim

Color : Clear

Heat Resistant

Chemical Resistant

Size : 16 x 150 mm

Thickness min: 1.0mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/