



Product Code . DS-LE-10919

Reagent Bottle with Glass Stopper

Description

DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass

Capacity : 100 ml, dia:52mm height: 97, Neck: 29/22

Capacity : 250 ml dia: 71mm height: 129 Neck: 34/24

Capacity : 500 ml dia: 86 mm height: 164 Neck: 45/40

Capacity : 1000 ml dia:107mm height: 200 Neck: 60/46

With glass stopper

Material : Borosilicate Glass

Heat resistant

Chemical resistant

Color : Clear

A maximum usage temperature of +100 °C is recommended.

Thermal shock resistance, 30 K. Hydrolytic class

If the bottle and the stopper have a temperature difference greater than 30 °C

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/