



Product Code . DS-LE-10924

## Electric Oven

### Description

---

#### DESCRIPTION :-

The Gap between the two walls filled with glass wool insulation.

Fully insulated door with brass lock hinges & aluminum ventilator.

Heating is by specially designed heaters made of Nichrome wire fitted at the bottom of the Oven.

Temperature Controlled by Electronic Digital Temperature controller with PT 100 as sensor.

Suitable for laboratory i/c one pair of heat resistant hand gloves

Digital Display

Capacity: Min 47 Liter

Temp. Range: Accuracy at Min  $120 \pm 0.3$

Uniformity Min:  $120 \pm 2.0$

Heat up time : Min 11



Shelves No.: Min 2 pcs,

Power Supply: AC110V/220V $\pm$ 10%,50/60Hz

---

Internal Size (W\*D\*H): Min 390\*290\*390mm

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)