

Product Code . DS-LE-10925

Viscometer Digital Display



Description

Description :-

Digital viscometer is a new digital product used for determining the liquid viscose capacity and the absolute viscosity.

Features:-

It is widely used to determine and measure the liquid viscosity in many applications such as grease, painting, pharmacy and adhesives.

High measuring accuracy

Stable in measured display

Easy operation and read-out

Excellent in Anti- interference.

Technical Specification : -

Measurement Range : 10-1X105 mPa.s

Rotor Speed (rpm) : 6/12/30/60

Rotor: 1, 2, 3, 4

Accuracy: $\pm 5\%$ Newtonian Fluid

Power: AC 220V $\pm 10\%$ 50 Hz $\pm 10\%$

Relative humidity no more than: 80%

Outline Dimension (mm): Min 98×135×165 mm

Net Weight (kg) Min 7 kg

Standard accessories with operation manual

Warranty for One Year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/