

Product Code . DS-LE-10926

Melting Point Apparatus



Description

DESCRIPTION :-

Easy-to-use, internal LED light

Digital display: the two different rise speeds allow for quick measurement.

Operating temperature °C: Min Room to +350

Precision at 20°C °C: ± 1.0

Precision at 300°C °C: ± 2.5

Display: Min 4 digit with 8 segments

Hold function: Yes

Number of samples: 3

Fast-speed ramp °C min: Min 20 fixed

Slow-speed ramp °C min: from 1 to 10

Capacity W: Max75,

External dimensions Length x Width x Height mm: Min 215x155x185

Weight: Min Kg 2

Standard accessories with printed operation manual.

Warranty for One Year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/