

Product Code . DS-LE-10927

Magnetic Stirrer & Hot Plate



Description

DESCRIPTION :-

Magnetic Stirrer With speed regulator, temperature regulator and mains, all stainless-steel body having stainless steel top magnetic stirrer with hot plates having different dimension.

Maintenance free brushless DC motor

Ceramic coated SS top plate

Aluminum casting body

Heating and stirring function

Speed: 100~2000 rpm

Durable stirring motor

Capacity: 5 liters

Power Consumption: 500 Watt



Max temperature: 380-degree Celsius

Plate dimension: 165 x 165 mm

Included with rod stand & magnet

Power source: AC220V±% 60Hz

Warranty for One Year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/