

Product Code . DS-LE-10928

Hand Held Refractometer



Description

DESCRIPTION :-

Measuring range: 0~32% sugar scale

With zero adjusting screw driver

Accuracy: $\pm 0.2\%$

EASY TO USE: Just take a drop or 2 of your brew on the prism and hold the refractometer to the light, then look through the eyepiece to get the reading.

No lookup table needed and no waste!

Dimension: min 25 x 38 x 139

High quality pure aluminum body construction

Supplied with box

Warranty for One Year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/