

Product Code . DS-LE-10929

Conductivity Meter



Description

DESCRIPTION :-

Designed in ABS Plastic to provide the rugged portability matched in fields with the accuracy demanded in the laboratory.

A beautiful carrying case holds everything needed in the field operation.

These are useful test instruments to measure Conductivity of any solution in the fields, laboratory or plant.

The portable conductivity meter provides fast, accurate readings with digital readability and the convenience of a remote probe 1 separately

Multi measuring ranges: 1.999 mS and 19.99 Ms

Resolution: 0.001 mS for 1.999 mS & amp; 0.01mS for 19.99 mS

Accuracy: $\pm (3\% \text{ F.S.} + 1 \text{ d})$

LCD display for low power consumption & amp; clear read-out even in, bright ambient light condition



Water resistance on the front panel

All function keys are used the rubber buttonUsed the durable, long-lasting components, including a strong, light

weight ABS-plastic housing case

Built-in Low battery indicator, Supplied with standard electrodes & accessories

Power supply: 9V battery

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/