



Product Code . DS-LE-10931

Atomic Model

Description

DESCRIPTION :-

A set of atomic models in plastic balls

These sets consist of moulded balls of different colors and sizes along with connecting lugs of different sizes.

These sets are very helpful in modeling different organic & Inorganic compounds.

Inorganic/organic molecular model student set includes: 51 atom-parts, 6 Carbon 4 holes tetra, black, 6 Oxygen, 2-holes ang., red, 14 Hydrogen, 1-hole, white, 3 Nitrogen, 2 tetra, 1 tribipyr., blue, 3 Sulphur, 1 ang., 1 tetra, 1 octa, yellow, 1 Phosphorus, 4-holes, purple, 6 Halogen, 1-hole green, 3 Metal, 1-hole, 17mm, grey, 6 Metal, 3 ang., 2 pyr., 1 tetra, grey, 3 Brown, 1 sp3, 1 dsp3, 1 d2sp3 26 links, grey

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/