



Product Code . DS-LE-10932

Platinum Wire

Description

DESCRIPTION :-

Mounted in rubber stopper

They are made in the form of a long glass rod with exposed platinum metal on the tip of the rod.

The rest of the conductor which is mainly inside the glass rod is made of a specific copper alloy.

Part of the conductor is also exposed for making electrical connection with your power source.

Diameter: 0.020 in. / 0.5 mm / 24 gauge

Overall length: 85 mm

Glass tubes, Warranty for One Year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/