

Product Code . DS-LE-10935

Thermometer



Description

DESCRIPTION :-

Measuring range: -40°C ~ 1200°C

Resolution: 1°C

Accuracy: $\pm (1\% + 1^\circ\text{C})$ for 0°C to 750°C

Display: 3 1/2 digits

Sensor type: Thermocouple K

Two channel input facility (T1, T2)

Sampling time: 0.4 second

Power source: DC 9 Volt battery

Supplied with temperature probe and printed operation manual

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/