

Product Code . DS-LE-10945

Burette Stand



Description

DESCRIPTION :-

A Retort Stand is a piece of scientific equipment to which clamps can be attached to hold other pieces of equipment and glassware for instance, burettes, test tubes and flasks.

Made of cast iron powder coated mild steel rod chrome plated.

Features:

Best quality raw material

Long life

Most appreciated

Made of best quality & heat proof Borosilicate glass



To hold one burette properly

Made of cast iron

Length: min 50 cm

Width: 60 cm

Heavy base & clamp

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/