

Product Code . DS-LE-10946

Thistle Funnel



Description

DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass

Material: Borosilicate glass.

Heat resistant

Chemical resistant

Outer Diameter of top : 38 mm

Outer Diameter of stem : 7 mm

Total Length : 350 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/