



Didac Scientific



Product Code . DS-LE-10947

Wall Room Thermometer



Description

DESCRIPTION :-

Made of fine quality plastic body

Temperature Range : Min- 580°F~ 158°F (-50°C -70°C)

Temperature Accuracy : ± 1 °F

Humidity Accuracy: 30%-80%

out of range $\pm 3\%$, All time high and low records



Refresh Rate : Min 10 Sec., °F/ °C Selector

Included Components : Temperature and Humidity monitor, AAA battery user manual. Distinctive moulding in a very large clear print that is easily read.

Graduated

Arrangement for hanging on wall

Warranty for One year.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/