



Product Code . DS-LE-10948

Desicator

Description

DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass

With porcelain plate

Capacity : min 180 mm

Operating Temperature: -25 to 120 Degree Celsius.

Depth to plate support : 100 mm

Total Height : 280 mm

Diameter of porcelain plate : 150 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/