

Product Code . DS-LE-10953

## Separating Funnel



### Description

---

#### DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass

---

Squib pear shape with ground in glass stopper and stopcock

Capacity: 250 ml



Height: 280 mm

Heat resistant

Chemical resistant

Color: Clear

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)