



Product Code . DS-LE-10964

Lab Table

Description

DESCRIPTION :-



Size : 2400(Length) x 1500(Width) x 900 mm(Height)

Made of scratchproof pesticide treaded melamine coated laminated board,

Table top 30 mm and others panels are of 18 mm thickness, Edges of panels and top are to be sealed 2mm thick PVC edging by automatic edging banding machine top and panels would be joined by using housing dowel, bolt, screw T-nut and pneumatic nailing where necessary and 4x2 drawer and 2 open shelf in both side and open both side open rack on the table top door for electric duct, switch board with cable passing system.

As per attached drawing, Detail Drawing would be submitted showing dimensions of the sections.

Color: Gray & Graphite

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/