

Product Code . DS-LE-10968

Cork Borer Set



Description

DESCRIPTION :-

Made of brass



Chrome plated finish with individual handles

Cutting edge at bottom with rods for cleaning borers

Each set of 6 nos

The sets come complete with push rod for clearing borers.

Diameter : 12 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/