



Product Code . DS-LE-10975

Bottle (Aspirator)

Description

DESCRIPTION :-

High Quality Aspirator Bottle made from Glass

With Interchangeable Stopper & Stopcock with Outlet for Tubing

Heat Resistant

Chemical Resistant

Color : Clear

Capacity: 5000 ml

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/