



Product Code . DS-LE-10976

Lighting Condenser

Description

DESCRIPTION :-

Made of best quality & heat proof Borosilicate glass With fused inner tube

Length of jacket : 400 mm

O.D. of upper tube : 20 mm

Material : Glass

Used in laboratory

Heat resistance

Extended durability

Total Length : 580 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/