



Product Code . DS-LE-10980

Glass Rods

Description

DESCRIPTION :-

Hard Glass and Best quality

A vast assortment of Glass Rod or Stirrers that is used for mixing chemicals when performing experiments.

Made using borosilicate glass, these glass rods are durable in nature.

These rods find wide application in laboratories.

Heat resistant

Chemical resistant

Color : Clear

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/