



Product Code . DS-B-10988

Glass Rod

Description

Catalogue No. : "IDAC5108/1.2"

DESCRIPTION :-

Material : Borosilicate Glass

Heat resistant

Chemical resistant

Color : Clear

Used for stirring solutions during experiments.

Size: Diameter X Length : 6 X 300mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/