

Product Code . DS-B-10990

Compound Microscope





Description

Catalogue No. : "IDAC7111/1.4"

DESCRIPTION :-

- Student Microscope
- Illumination by plano concave mirror in fork mount.
- Steple Grey, chemical resistant, back-on finish.
- U-shaped, well balanced base with heavy inclinable body. It has two spring clips, built-in condenser, disc-diaphragm.
- It has provision for fitting mechanical stage.

-
- Separate coarse and fine focusing is provided by rack and pinion arrangement.
 - (Triple revolving nosepiece has accurate click stop arrangement for perfect centering and identical repeated settings).
 - Objectives 5x,10x and 45x, Eyepieces 10x, 15x
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/