



Product Code . DS-B-10996

Blade Scalpel Handle

Description

Catalogue No. : “IDAC4223/1.10”

DESCRIPTION :-

- Scalpel Handle
- Material : Stainless Steel
- Usage : Hospital, Clinical
- Type : Non-Disposable
- It is made of Stainless-Steel having size no. 4.
- It can be used with blade with blade no. 20, 21, 22, 23, 24, 25. (Only handle without blade)

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/