



**Product Code . DS-B-11002**

## Beaker

### Description

---

**Catalogue No. : “IDAC1140/1.16”**

**DESCRIPTION :-**

- Material : Polypropylene
- Capacity : 250ml
- Beakers, moulded in Polypropylene
- Excellent clarity and very good chemical resistance.
- These autoclavable Beakers have prominently raised graduations on the outer surface.
- Pouring a solution becomes very convenient because of its tapered spout.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)