



**Product Code . DS-B-11030**

## **Needle & Brush**

### **Description**

---

**Catalogue No. : “IDAC7474/3.12”**

#### **Dissecting Needle :-**

- Handle Material : Plastic
- Chemical resistant
- This is a straight needle with plastic handle for general dissection use.
- It's a versatile tool with many uses in the laboratory.

#### **Brush :-**

- Usage/Application : Cleaning
- Material : Nylon
- Bristle Stiffness : Soft
- Shape : Round
- Color : White
- Used for cleaning tubes.

