



Product Code . DS-B-11033

Beaker

Description

Catalogue No. : “IDAC1284/3.15”

DESCRIPTION :-

- Material : Polypropylene
- Capacity : 250ml
- Beakers, moulded in Polypropylene
- Excellent clarity and very good chemical resistance.
- These autoclavable Beakers have prominently raised graduations on the outer surface.
- Pouring a solution becomes very convenient because of its tapered spout.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/