



Product Code . DS-B-11034

Weighing Balance

Description

Catalogue No. : “IDAC2987/3.16”

DESCRIPTION :-

- Triple Beam Balance are specially designed to handle weighing solids and powder.
- Three notched and tiered beams with center reading.
- Indicating sliding masses giving total weighing capacity of 610 gm.
- The capacity can be increased up to 2610 gm. by using supplementary masses.
- Spring loaded zero adjust compensator and magnetic damping are provided for fast setting and reading..
- Capacity : 610 gm.
- Readability : 0.1 gm.
- Front Beam Calibration : 10 x 0.1 gm.
- Second Beam Calibration : 500 x 100 gm.
- Rear Beam Calibration : 100 x 10 gm.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/