



**Product Code . DS-B-11042**

## Blade

### Description

---

**Catalogue No. : “IDAC1211/5.2”**

**DESCRIPTION :-**

- Scalpel Handle
- Material : Stainless Steel
- Usage : Hospital, Clinical
- Type : Non-Disposable
- It is made of Stainless-Steel having size no. 4.
- It can be used with blade with blade no. 20, 21, 22, 23, 24, 25. (Only handle without blade)

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)