




Product Code . DS-B-11045

Wire Gauze

Description

Catalogue No. : "IDAC4522/5.5"

DESCRIPTION :-

- Size : 15 X 15cm
- Flat gauze ceramic filled centre for use on tripods.
- Wire Gauge With asbestos center, extra strong quality in squares.
- Shape : Square

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/