



Product Code . DS-B-11047

Cover Slip

Description

Catalogue No. : "IDAC6102/5.7"

DESCRIPTION :-

- Perfect for preparing permanent mounts.
- Size 18x18mm
- Square Shape.
- These microscope slide cover slips are made out of glass and are essential equipment for the lab.
- They can be used to keep solid specimens pressed flat, and liquid samples shaped into a flat layer of even thickness.
- Material: Borosilicate Glass
- Color: Clear
- Used in Laboratory

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/