



**Didac Scientific**



**Product Code . DS-B-11050**

## **Measuring Cylinder**



---

## Description



---

**Catalogue No. : “IDAC9611/5.10”**

### **DESCRIPTION :-**

- Material : Polypropylene
- The Measuring Cylinders, moulded in Polypropylene & TPX, are single piece cylinders.
- These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.
- The broad base provides good stability.
- Capacity: 250ml.
- Least count: 2ml

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)