



Product Code . DS-B-11058

Hot Plate

Description

Catalogue No. : "IDAC0237/5.18"

DESCRIPTION :-

- Diameter: 200mm
- Shape: Round
- Electrically operated.
- Frequency (Hz): 50Hz
- Material: Stainless Steel
- Round Plate Top
- Plate Diameter 20 cm.
- With Energy Regulator and ON-OFF switch
- Working voltage is 230V AC 50Hz

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/