



Product Code . DS-B-11063

Bunsen Burner

Description

Catalogue No. : "IDAC9844/6.3"

DESCRIPTION :-

- Burner bunsen without stopcock
- Common piece of laboratory equipment that produces a single open gas flame, which is used for heating, sterilization, and combustion.
- Nickel plated burner tube of 12mm dia with air regulator and a riffled connector.
- Mounted on a painted cast iron base.
- Only the air supply port can be adjusted. without stopcock

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/

