

**Product Code . DS-B-11067**

## **Conical Flask**



### **Description**

---

**Catalogue No. : “IDAC9162/6.7”**

**DESCRIPTION :-**

- Capacity: 100ml
- Material: Polypropylene
- These cone shaped flasks are rigid, translucent & autoclavable.
- The screw cap keeps the contents of the flask contamination free and makes the flask leakproof.
- Threads disengaged before autoclaving.

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)