



Didac Scientific



Product Code . DS-B-11069

Measuring Cylinder





Description

Catalogue No. : “IDAC2174/6.9”

DESCRIPTION :-

- Material: polypropylene
- The Measuring Cylinders, moulded in Polypropylene & TPX, are single piece cylinders.
- These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.
- The broad base provides good stability.
- Capacity: 250ml.
- Least count: 2ml

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/