



Product Code . DS-B-11075

Test Tube Stand

Description

Catalogue No. : "IDAC3399/7.4"

DESCRIPTION :-

- Test Tube Stand Plastic 6 Holes 6 Pegs
- This autoclavable Stand is perhaps the most commonly used Test Tube Stand in various laboratories. For better visibility, this 6-place rack is end supported by 2 columns.
- The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes whereas the base has corresponding hemispherical wells.
- It also has six vertical pins for drying tubes. Another option available in this section is a test tube stand in which all the holes in the top plates can hold 25 mm dia tubes.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/