

Product Code . DS-P-11110

Universal Clip



Description

Catalogue No. : “IDAC0555/1.6”

DESCRIPTION :-

-
- Made of Aluminum-alloy die casting.
 - This Laboratory clamp has PVC coated jaws capable of holding upto 90 mm article.
 - It is provided with 8mm diameter & 6" long Stain less Steel rod.
 - It is coated with LEAD-FREE powder coating & fitted with brass screw.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/