

**Product Code . DS-P-11123**

## **Potentiometer 4 Wire**



---

## Description



---

**Catalogue No. : “IDAC7151/2.1”**

**DESCRIPTION :-**

- Potentiometer 4 wire with pencil jockey.
- Four constantan wires 24 swg stretched on a polished wooden base connected in series with the help of brass strips provided with two 4mm socket terminals.
- A box wood meter scale fitted on one edge of the wooden base parallel to one meter length of 4 wires with a knife edge jockey.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)