

Product Code . DS-P-11127

Galvanometer



Description

Catalogue No. : "IDAC7665/2.5"

DESCRIPTION :-

- Centre Zero Universal Meter Galvanometer
- The centre zero universal meter (Universal Galvanometer) is electrically identical to the universal meter.
- 1000-ohm movement, moving coil, diode protection against overload up to 12V D.C.
- Range: DC 30-0-30 (20 μ A/ DIV)
- The unit is ideal for use as an indicator in slow A.C.
- Housed in ABS plastic molded case

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/