



Product Code . DS-P-11132

Ammeter

Description

Catalogue No. : "IDAC3589/2.10"

DESCRIPTION :-

- Range DC: 0-10 amp
- An Ammeter is a measuring instrument used to measure the current in a circuit.
- Electric currents are measured in amperes, hence the name.
- Instruments used to measure smaller currents, in the milliampere or microampere range, are designated as milliammeters

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/