



**Product Code . DS-P-11141**

## **Meter Bridge**

### **Description**

---

**Catalogue No. : “IDAC8234/3.1”**

#### **DESCRIPTION :-**

- The four-gap bridge has large and heavy brass strip mounted on a heavy  $\frac{3}{4}$  thick laminated base of size 44 x 5 complete with pencil jockey and lead for connection.
- Dimension: 1.10 m x 100 mm x 50 mm high including terminals

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)